

StackShip

IATA Approved Transport Shippers

- Fully autoclavable multiple times
- Double wall polypropylene insulates totally stabilises temperatures
- Light weight material helps keep freight cost down

OpenD

Your Scandinavian Distributor
www.opend.eu



A wire lined viewing window ensures safe viewing and no escapees

IATA approved markings with required warning text

Double mesh protected filter panel to ensure filters are never damaged or compromised



Up to 5 dividers fitted into either size shipper




Tapered sidewalls means filter area can't be blocked and boxes are stackable. Delivered assembled to save time.

RECYCLING POSSIBILITIES:


Recycle your shippers into Surgery Pads. Ideal for all work including necropsy's etc. Simply cut 2 pads between 10 and 20cm square, tape the 2 layers together and you have a very durable surgery pad that won't absorb blood like cork pads, is very easy to clean and wipe down, is autoclavable and lasts a long time.

Use our protective absorbant PuraMat with this surgery surface.

You can get between 2 and 4 surgery pads from every shipper!




cut 2 pads between 10 and 20cm square



tape the 2 layers together



durable surgery pad, easy to clean, accept pins, autoclavable



Delivered stackable and assembled to save you time and cost

- 170 small / pallet
- 100 large / pallet



Quick strap for greater safety

Stackable Shippers

- IATA Transport Authority Approved
- Polypropylene Shippers Autoclavable @ 121°C
- Rat, mouse, guinea pig and all small animals
- Delivered Assembled with Internal Viewing Windows



Large Animal Transport Shipper 580 x 370 x 210mm H (Floor Area: 1300cm²)

Capacity: 30 mice up to 25gm. 5 rats up to 250gm maximum Internal Height: 200mm

Code	Description	Pack of	Price/Fr	
			20x	or 100x
AS1300ASTSAS	Assembled Shipper with Lid & Window	100	411	1,710
AS1300LD	Loose Dividers	100	115	480
AS1300FCS	Fly Wire Cover Seal (For Lid)	100	125	520
ASVW	International View Window	100	-	
ASNECTAR	Necta Moisture Sachets 100 x 75g	1	154.00/box	



Small Animal Transport Shipper 410 x 320 x 150mm H (Floor Area: 750cm²)

Capacity: 12-15 mice maximum Internal Height: 150mm

Code	Description	Pack of	Price/Fr	
			20x	or 150x
AS750ATSAS	Assembled Shipper with Lid & Window	150	369	2,303
AS750LD	Loose Dividers	150	108	675
AS750FCS	Fly Wire Cover Seal (For Lid)	150	108	675
ASVW	International View Window	150	-	
ASNECTAR	Necta Moisture Sachets 100 x 75g	1	154.00/box	



Rabbit & Poultry Transport Shipper

- IATA Approved
- Polypropylene Coreflute Repeatedly Autoclavable @ 121°C
- Poultry, Ferret, Rabbit and Cats
- Wire Lined Filtered Windows
- Floor Area: 2450cm²
- Dimension: 600 x 410 x 310 H mm
- Internal Height 300mm

Code	Description	Packs of	Price/Fr	
			20x	or 120x
AS2200PTS	Poultry/rabbit Transport Shipper Flat Pack	10	810	4,050
	Bulk Pallet Qty's	120	-	



Don't forget your Necta H20 Gel

Staplers, Hand Strappers and Clips as well as Automated Strapping machines for bulk shipping available. Please inquire.

