

## MOJO-CLIP™

Making life easy, economical and simple. **Get your Mojo.**

- Makes sample taking, Biopsy and DNA simple and fail-safe
- The 'Easy-Clip' thumb punch mechanism has been incorporated into the design of the MojoClip™ to allow for easier use and prevention of RSI's (Repetitive Strain Injuries).
- The MojoClip™ can be used as a method of animal identification or for the collection of DNA samples with Eppendorf Caps.
- Each MojoClip™ is individually sharpened and checked by the Onyx Quality Assurance department before sale, giving you peace of mind that you are using another quality ONYX instrument. 100% Satisfaction Guaranteed. \*
- Low-Cost Automated Production, All stainless steel, non-corrosive




**Code:** 251-1094

**Description:** Premium Onyx 'easy-clip' punch 7cm, 2mm hole, with Eppendorf holder

**Purchase Quantities:**

1-5	
5-10	
10-50	
50+	Request Offer



**Key words:**

Ear Punch, Tail, Foot, Biopsy, DNA, Eppendorf, Sampler, Clipper, Sharpened, Warranted, Simple, Quick, Easy, Low Cost, Quality,