

Surgical Instruments

The ear punchers are all autoclavable but should be lubricated afterwards for a long and trouble free usage.



MOJO-CLIP™

Code 251-1094

Description Premium Onyx 'easy-clip' punch 7cm, 2mm hole, with Eppendorf holder



Premium ONYX Thumb Punch

Code S00076

Description Premium ONYX notcher with unique 'easy-clip' thumb lever. 7 cm, 2 mm hole. Individually Sharpened, Checked and Warranted.



S/S Thumb Punch

Code S500075

Description S/S Thumb Punch. Individually sharpened, checked and warranted.



Ear punch pincer

Code AT7000, hole 2mm OD Stainless Steel

Code AT7071, hole 1mm OD Stainless Steel on request

Code AT7075, hole 2,75mm OD Stainless Steel

Description After hearing several testimonies from customers with sore wrists and thumbs from using small ID punchers, we developed this ear punch pincer. The AT7000 pincer relieves the stress on the wrist during ID marking as it requires very gentle pressure to cut a hole. The stainless steel puncher has a sloped edge which provides a cut rather than a punch, which minimizes animal discomfort.



Ear punch pincer

Code AT7020, hole 2mm OD Stainless Steel.

Description Small and efficient ear puncher based on the AT7000. The puncher allows the user to come closer when marking the animal, and it is also equipped with a sloped edge for a clean cut.