

*Making transfer
of cage from rack
to trolley easy!*



Height Adjustable Trolleys

– For transport and lifting of heavy and bulky items or for use as work station

These trolleys offer assistance with lifting and transport of heavy and bulky items:

- Ideal for lifting and transport of large cages, bottle baskets, etc.
- Large trolley is also ideal for procedure work and cage change ...
- Available in two standard sizes

Large trolley: (shown on brochure)

Work surface dim: 1250 x 525mm

Height adjustable from 720-1390mm

• Load: 140 kg

Small trolley: (not shown - on request)


Work surface dim: 625 x 525mm

Height adjustable from 853-1353mm

• Load: 70 kg

OpEnd

Life Science Solutions



Work platform is 100% flush with the surface, no edges or bumps.

Easy to adjust the height manually or use the programmable height settings.

Three programmable standard heights e.g.:

Setpoint 1—activate this button:

The trolley platform will stop at the height programmed into this button e.g. top shelf of rack.

Setpoint 2—activate this button:

The trolley platform will stop at the height programmed into this button e.g. your preferred working height.

Optional features:

- Battery operated
- Container / Pencil holder with perforated bottom
- Drawer with perforated bottom
- Lower shelf

OpEnd

Life Science Solutions

OpEnd
DK-4681 Herfølge, Denmark
T: +45 29 344 334
E: tine@opend.eu
W: www.opend.eu