



Customized to your specific requirements:

- Cage change
- Cage cleaning out
- Dosing
- Surgery
- Etc...



OPendAIR
– Clean Air Solutions

Laminar Air Flow Cabins

– For Lab Animal Science Work

Specially designed to meet the requirements in LAS

Ensuring a clean HEPA filtered working environment with a one way laminar air flow providing protection of:


- the operator
- the environment
- the animals

Various sizes and options are available:

- Automatic cage dispensers
- Height adjustable working table
- Safe work with anesthesia solution
- Cage cleaning out solution
- Fixed / Sliding curtains
- Solid curtains / strip curtains
- Outlets for power, data, gas, etc
- PC / Tablet holder

The LAF cabin is tested in accordance with the following directives:

- ISO 5, jf. DS/EN ISO 14644-1 Classification of air cleanliness by particle concentration (Cleanrooms and associated controlled environments)
- EN12469 Biotechnology – performance criteria for microbiological safety cabinets.
- Sterile Drug Products Produced by Aseptic Processing 2004



Heated work surface in Corian



Laptop holder



Designed for safe work with anesthesia



Red light



White light

We make customized solutions – meeting your specific requirements!